

wheelchair stickers

Fitting your wheelchair stickers

A simple guide for your spoke guards

Before you start: use a clean spray bottle filled with water and a couple of drops of washing-up liquid.

Stickers + guards? We fit your stickers before posting.



1 Mist the guard

Lightly spray the surface of your plastic spoke guard.



2 Peel & mist

Leave the rectangular backing sheet flat. Peel up one of the two stickers, then lightly spray its grey adhesive side.



3 Position your design

Place the sticker onto the guard and slide it into position.



4 Smooth it down

Smooth the sticker with a clean towel to remove excess water.



Fasten the guards to your wheels

Three close pairs of small holes. Feed each of the supplied clear cable ties from the back through one hole and back through its neighbour. Fasten around the nearest spoke behind the guard, keeping the locking head at the back.

The small fitting cut lets the sticker overlap slightly over the guard's domed shape.

www.wheelchairstickers.com | info@wheelchairstickers.com

Illustrations show the fitting method. Your artwork will vary.